

# Phonak CROS P Order Form

Phonak CROS P is a solution for unilateral hearing loss designed to help the wearer follow conversations. CROS P is worn on the unaidable ear to pick up sounds and wirelessly transmits them to a Phonak Audéo™ P hearing aid worn on the ear with better hearing. CROS P has no acoustical output.

CROS P includes two models and is compatible with Audéo P-13T, Audéo P-R and Audéo P-RT hearing aids in all performance levels.

## Step 1: Order details

Hearing Clinic:	
Delivery Address:	
Contact Name:	Phone Number:
Email Address:	
Date Ordered:	Date Required:

## Step 2: Patient information

Name: \_\_\_\_\_

## Step 3: Choose hearing instrument

Model	Part No.	Colour	Qty
CROS P-R	050-0860-xx		
CROS P 13	050-0861-xx		


## Colour chart



## Step 4: Choose your Acoustic Coupling Option

### CROS Slim Tube with Domes


CROS Tube is available in left/right and four different lengths (0, 1, 2, 3) with standard dome available in three sizes (S, M, L)



L R Select options in section a and b only below

### CROS Slim Tube with Slim Tip\*

CROS Tube is available in left/right and four different lengths (0, 1, 2, 3), with a custom made earpiece



L R Select options in section b and d only below

^ Impression required

### CROS Tip\*

Fully custom made earpiece for CROS, available in four different cable lengths (0, 1, 2, 3)



L R Select options in section c and d only below

^ Impression required

## Step 5: Choose from the following options

### (a) Domes 4.0

Dome style	Dome size	
	Left	Right
Cap	Left	Right
Open	Left	Right
Vented	Left	Right
Power	Left	Right

### (b) CROS Slim Tubes



### (c) CROS Tip Cables - must be paired with a CROS Tip below

054-0932	CROS Tip Cable 4.0 0L	054-0933	CROS Tip Cable 4.0 0R
054-0934	CROS Tip Cable 4.0 1L	054-0935	CROS Tip Cable 4.0 1R
054-0936	CROS Tip Cable 4.0 2L	054-0937	CROS Tip Cable 4.0 2R
054-0938	CROS Tip Cable 4.0 3L	054-0939	CROS Tip Cable 4.0 3R

### (d) Slim Tip 4.0 or CROS Tip



### (e) Options

- SlimTip Acrylic (062-0006-01)
- SlimTip Silicone (062-0007-01)
- CROS Tip 4.0 Hard Acrylic (062-0016-\*\*) - must be paired with a CROS Tip cable

\*Open Vent for SlimTip by default

### (f) Generic options (all earpieces)

- Removal line
- Removal line extension + 5mm
- Locks